

---

Subject: Re: prespec data analysis package  
Posted by [miree](#) on Wed, 23 Jul 2014 08:37:28 GMT  
[View Forum Message](#) <> [Reply to Message](#)

---

Hi,  
I've never seen this type of background around the peak in a MHTDC spectrum... I cannot say what is causing it.

The timestamp difference histogram looks good, except that I've never seen this structure close to 0. Did you have the target in the forward (high intensity) position?

Perhaps you can make a trigger selection when looking at the spectra. How does the histograms look with a particle trigger?

You can add the line  
triggers 9 10  
inside any processor and it will only process data and display histograms for events with trigger 9 or trigger 10.

---