
Subject: Re: TG4SDServices::GetVolumeID:
Posted by [asanchez](#) on Mon, 04 Jun 2007 09:34:37 GMT
[View Forum Message](#) <> [Reply to Message](#)

Hi

I have removed the underscore symbol from the sensitive detector name, and it seems to work.

But when I run now, the PndBoxGenerator for a multiplicity of 100 particle(whatever..) i get this error information:

>>> End of Event 0

-E- CbmStack: Particle index -999 not found in dex map!

Fatal in <CbmStack::CbmStack::UpdateTrackIndex>: Particle index not found in map aborting

Using host libthread_db library "/lib/libthread_db.so.1".

Attaching to program: /proc/28321/exe, process 28321

`system-supplied DSO at 0xffffe000' has disappeared; keeping its symbols.

done.

done.

[Thread debugging using libthread_db enabled]

[New Thread 16384 (LWP 28321)]

Any ideas?

Thanks a lot
alicia.
