

---

Subject: Re: problem with TPC recotasks

Posted by [Anonymous Poster](#) on Thu, 27 Mar 2008 13:01:32 GMT

[View Forum Message](#) <> [Reply to Message](#)

---

Hi Alicia,

Just changing (id>100000) to (id>1) doesn't work. If I do that, then after reconstruction most of the events have only one track or even Zero. Could you please let me know how to get the only the primary particles (for Lambda0 case only proton and pion).

What does the "unsigned int id=point->GetTrackID();" do? In some cases this id are >2000 or >3000. What are these values?

Best Regards,

Dipak

---